

PRODUCT CATALOGUE

2020-2021



Compression, Vascular, Lymphoedema,
Wound Care Products.



LEADER IN PHLEBOLOGY



A LONG EUROPEAN HISTORY

We are proud to present our products for the compression therapy; they truly represent the essence of our company and stem from a tradition handed over from father to son. This century-long involvement has lead us, year after year, to develop medical stockings that match advancements in modern phlebology with state of the art technology, offering products that provide the optimum medical effect. Gloria's history began in 1895 as medical business in Milan, Italy, while manufacturing activity started in 1926 in Lugano, Switzerland. In 1936 appeared the first, totally automated loom to produce two-way stretch stockings. In the 50's Gloria moved to Menaggio, on lake Como becoming a well-known brand in the phlebology field. In 1990, Gloria joined the famous Quality Association of medical stockings producers in Germany. Today, it is among the leaders in its field, with the quality of its products recognized all over the world.

C O N T A C T S



(+44) 01462 411819



info@gardamed.com



www.gardamed.com



GLORIA MED GROUP®

Gloria Med S.p.A.: Via Diaz 7 - 22017 Menaggio (CO) - ITALY

GLORIA MED

6

Production process	6
Quality certifications	6
Gloria compression therapy	7
Continuos evolution	7

COMPRESSION THERAPY

8

Medical Compression stockings

Gloria Med Gentleman	8
Gloria Med Micro	10
Gloria Med Soft	12
Gloria Med Comfort	14
Gloria Med Cotton	16
Gloria Med Strong	18
Gloria Med Perfect	20
EasyVen	23

Venous ulcer

UlcerKit	24
UlcerKit PRO	26

Anti-embolism stockings

Gloria Med Thrombosis	28
Post Op	29

Compression Arm Sleeve

Gloria Med Sleeve	30
Gloria Med Perfect	32

Support Stockings

Lady Gloria	34
Man'Support	38
Man'Support Silver Cotton	40
Diabesoft	41
Travel & Work	42

DERMA COSMETICS PRODUCTS FOR HEALTHY LEGS

43

HappyVen Fresh/Soft	43
---------------------------	----

PRODUCTS FOR PHYSICIAN

44

Sclerotherapy

Siderlink	44
-----------------	----



PRODUCTION PROCESS

The effectiveness of the compression therapy, and thus the well-being of the patient, is a consequence of top quality products; Gloria stockings are manufactured with the best yarns, careful processing, computer-controlled knitting-machines, precise finishing and undergo final inspections. The results of this demanding process are two-way stretch, seamless stockings, perfectly matching the demand of modern phlebology while continuously undergoing extensive refinements in our R&D division .



QUALITY CERTIFICATIONS



Gloria Stockings are Class 1 Medical devices and meet all requirements of the EU Directive 93/42/EWG.



Gloria Compression stockings meet the RAL-GZG 387 requirements for German Quality Seal.



Stockings are routinely tested by the Hohenstein Research Institute in Germany



Gloria internal quality control system is UNI EN ISO 9001 and ISO 13485 certified.



Gloria stockings are FDA approved




Gloria Med is part of Eurocom, the "Association of European Manufacturers of Medical stockings with RAL certification and orthopaedic braces.

GLORIA COMPRESSION THERAPY

Compression therapy is used in the treatment of venous diseases and lymphedema. The aim of treatment is to maintain a high quality way of life despite an insufficient venous system. The type and extent of treatment depends on the severity of the problem and a doctor should always be consulted.

Gloria medical stockings are prescribed by physicians throughout the world, come in a wide selection of fabrics, styles, sizes, colours and compression ranges.



HOW DOES COMPRESSION WORK?

Compression therapy works by applying controlled pressure to the surface veins, keeping their diameter small and forcing blood back into the deep vein system which in turn pushes the blood uphill towards the heart. Stronger compression is applied at the ankle and decreasing up the leg. Improved blood flow aids in the healing of varicose veins, venous leg ulcer and lymphatic problems.

THE BENEFIT OF GLORIA COMPRESSION THERAPY

- Great for tired, achy legs
- Helps reduce and prevent swelling
- Assists with venous return of blood flow
- Recommended for venous insufficiency and lymphedema
- Maintaining a good quality way of life

CONTINUOUS EVOLUTION

Continuous evolution of scientific research, progress in technology and ever-growing demands by patients dictate that a modern industrial group must continually strive for improvements. Participation to international congresses and trade fairs, constant ties with doctors, research institutions, orthopaedical engineers and medical experts are a critical factor to our success. This is what drives Gloria and why its products are appreciated all over the world.





MICROFIBER RIBBED DESIGN

CL. 1 & 2

EVERYDAY COMFORT AND ELEGANCE ALSO FOR MEN

Medical compression socks made of thin microfiber yarn that combined with a modern ribbed design, highlights the male appearance and makes the product soft and elegant without losing the therapeutic effects of compression.

Available in three stylish colors that make these classic men's socks extremely versatile and comfortable.

COMPRESSION CLASSES:

- Gentleman 199: compression class 1
- Gentleman 299: compression class 2


INDICATIONS:

Class 1: (18-21 mmHg)


Mild symptoms of venous diseases such as heaviness and fatigue in the legs, mild varicosities without significant edema, post sclerotherapy, prophylaxis of thrombosis.

Class 2: (23-32 mmHg)

Moderate chronic venous insufficiency (CVI), thrombophlebitis (superficial thrombosis) with or without edema, to guarantee the therapy success (post sclerotherapy treatment, ambulatory selective varicose vein ablation, stripping) prophylaxis of thrombosis and embolisms.



AVAILABLE STYLES WITH CLOSED TOE



AD

KL.1 199


KL.2 299

(60) Black

(62) Dark Grey

(72) Blue Marine

MEASURING POINTS



SIZES (cm)

S - M - L - XL
NORMAL - MAX
SHORT - LONG

CIRCUMFERENCE MEASUREMENTS (cm)

	S		M		L		XL	
	normal	max	normal	max	normal	max	normal	max
d	31-34	35-40	35-38	39-44	39-42	43-50	43-47	48-55
b	20-22		23-25		26-28		29-31	

LENGTH MEASUREMENTS (cm)

	SHORT	LONG
A-D	34-38	39-44

COMPOSITION: Polyamide - Elastane

RAL 387 certified
medical stockings





MICROFIBER CL. 1 & 2

SHEERNESS, ELEGANCE AND MEDICAL COMPRESSION

Stockings and pantyhose suitable for women who look for softness and elegant appearance, with medically effective compression.

The use of state-of-the-art technologies and highly elastic yarns guarantee the best possible comfort and fit.

COMPRESSION CLASSES:

- Micro 181/191: compression class 1
- Micro 281/291: compression class 2

INDICATIONS:

Class 1: (18-21 mmHg)

Prophylaxis of venous diseases, heaviness and fatigue in the legs, pain with or without varicosity, initial varices during pregnancy.


Class 2: (23-32 mmHg)

Moderate CVI, post-varicophlebitis, venous insufficiency due to varicosis, after small superficial ulcer, post-sclerotherapy, varicosis during pregnancy.



AVAILABLE STYLES WITH
CLOSED TOE

OPEN TOE




AD



AGH*



AT



AD



AGH*



AGG**



AT

	AD	AGH*	AT	AD	AGH*	AGG**	AT
KL.1	191	191	191	181	181	181	181
KL.2	291	291	291	281	281	281	281

(34) Skin

(58) Mink

(60) Black

(72) Blue Marine


(62) Dark Grey

(68) Safari

* AGH with stay-up silicon band

** AGG: adjustable waist belt, left or right style.

MEASURING POINTS



SIZES (cm)


XS - S - M - L - XL
NORMAL - MAX
SHORT - LONG

CIRCUMFERENCE MEASUREMENTS (cm)

	XS		S		M		L		XL	
	normal	max	normal	max	normal	max	normal	max	normal	max
g	42-46	47-54	47-51	52-60	52-57	58-67	58-63	64-74	64-70	71-83
d	28-30	31-36	31-34	35-40	35-38	39-44	39-42	43-50	43-47	48-55
b	18-19		20-22		23-25		26-28		29-31	

LENGTH MEASUREMENTS (cm)

	SHORT	LONG
A-D	34-38	39-44
A-G	62-71	72-82



Lace stay-up band
Standard



Stay-up band with
silicon drops
On demand



COMPOSITION: 60% Polyamide - 40% Elastane



SYNTHETIC CL. 1 & 2

GRADUATED MEDICAL COMPRESSION WITH ELEGANCE

Stockings and pantyhose where compression is combined with a delightful appearance. Fine yarn and availability of open or closed toe to provide optimum comfort and satisfaction for patients. Gloriamed Soft are particularly suitable for young and dynamic women who need a therapeutical effect without compromising elegance and style.

COMPRESSION CLASSES:

- **Soft 141/151:** compression class 1
- **Soft 241/251:** compression class 2

INDICATIONS:

Class 1: (18-21 mmHg)

Mild symptoms of venous diseases such as heaviness and fatigue in the legs, combined venous and arterial insufficiency, mild varicosities without significant edema, post sclerotherapy, small varicose dilations of cutaneous veins, initial varices during pregnancy.


Class 2: (23-32 mmHg)

Moderate symptoms of venous diseases such as aching heaviness and fatigue in the legs, chronic venous insufficiency, post sclerotherapy, prophylaxis of thrombosis, moderate to severe varicosities, post phlebitic syndrome, lymphedema, stasis dermatitis due to venous insufficiency, post fracture or post surgical edema, edema from burns, varices during pregnancy, pregnant patients with previous phlebitis.



AVAILABLE STYLES WITH
CLOSED TOE

OPEN TOE




AD



AT



ATM



AD



AG*



AGG**


	AD	AT	ATM	AD	AG*	AGG**
KL.1	151	151	151	141	141	141
KL.2	251	251	251	241	241	241

(34) Skin

(60) Black

* AG: available also with silicon band
** AGG: adjustable waist belt, left or right style.

MEASURING POINTS



SIZES (cm)

XS - S - M - L - XL
NORMAL - MAX
SHORT - LONG

CIRCUMFERENCE MEASUREMENTS (cm)

	XS		S		M		L		XL	
	normal	max	normal	max	normal	max	normal	max	normal	max
g	42-46	47-54	47-51	52-60	52-57	58-67	58-63	64-74	64-70	71-83
d	28-30	31-36	31-34	35-40	35-38	39-44	39-42	43-50	43-47	48-55
b	18-19		20-22		23-25		26-28		29-31	

LENGTH MEASUREMENTS (cm)

	SHORT	LONG
A-D	34-38	39-44
A-G	62-71	72-82



Stay-up band with silicon drops

RAL 387 certified
medical stockings



COMPOSITION: 62% Polyamide - 38% Elastane





FINE COTTON CL. 1 & 2

FINE COTTON COMFORT AND ELEGANCE

A complete range of stockings where superior quality cotton, excellent smoothness and comfort are combined with a delightful look and therapeutical effect. Hypoallergenic and transpiratory, they are well suited for those who suffer from dermatological problems. They are also ideal for warm and humid environments and for people with allergies to synthetic fibers.

COMPRESSION CLASSES:

- Comfort 142/152: compression class 1
- Comfort 242/252: compression class 2

INDICATIONS:

Class 1: (18-21 mmHg)

Mild symptoms of venous diseases such as heaviness and fatigue in the legs, combined venous and arterial insufficiency, mild varicosities without significant edema, post sclerotherapy, small varicose dilations of cutaneous veins, initial varices during pregnancy.

Class 2: (23-32 mmHg)


Moderate symptoms of venous diseases such as aching heaviness and fatigue in the legs, chronic venous insufficiency, post sclerotherapy, prophylaxis of thrombosis, moderate to severe varicosities, post phlebitic syndrome, lymphedema, stasis dermatitis due to venous insufficiency, post fracture or post surgical edema, edema from burns, varices during pregnancy, pregnant patients with previous phlebitis.



AVAILABLE STYLES WITH

CLOSED TOE


OPEN TOE



AD



AGH



AD



AF*



AG*



AGG**



AT



ATM

KL.1	152	152	142		142	142			
KL.2	252	252	242	242	242	242	242	242	242

(34) Skin


(60) Black

* AF style available only in class 2 normal or with stay-up silicon band.

* AG style: class 1 only with stay-up silicon band; class 2 normal or stay-up silicon band

** AGG: adjustable waist belt, left or right style.

MEASURING POINTS



SIZES (cm)


XS - S - M - L - XL
NORMAL - MAX
SHORT - LONG

CIRCUMFERENCE MEASUREMENTS (cm)

	XS		S		M		L		XL	
	normal	max	normal	max	normal	max	normal	max	normal	max
g	42-46	47-54	47-51	52-60	52-57	58-67	58-63	64-74	64-70	71-83
f	37-40	41-47	41-45	46-52	46-50	51-59	51-56	57-65	57-62	63-73
d	28-30	31-36	31-34	35-40	35-38	39-44	39-42	43-50	43-47	48-55
b	18-19		20-22		23-25		26-28		29-31	

LENGTH MEASUREMENTS (cm)


	SHORT	LONG
A-D	34-38	39-44
A-F	52-59	60-68
A-G	62-71	72-82



Stay-up band with silicon drops



RAL 387 certified medical stockings



COMPOSITION: Comfort 142/152: 52% Polyamide - 26% Elastane - 22% Cotton
Comfort 242/252: 44% Polyamide - 30% Elastane - 26% Cotton





COTTON CL. 1 & 2

100 % NATURAL COTTON FIBERS GENTLY CONTACT THE SKIN

Stockings and pantyhose with high cotton content for a comfortable contact with the skin. Hypoallergenic and transpiratory, the stockings are well suited for those who suffer from dermatological problems. They are also ideal for warm and humid environments, for elderly people, and for people with allergies to synthetic fibers.

COMPRESSION CLASSES:

- **Cotton 162:** compression class 1
- **Cotton 262:** compression class 2

INDICATIONS:

Class 1: (18-21 mmHg)

Mild symptoms of venous diseases such as heaviness and fatigue in the legs, combined venous and arterial insufficiency, mild varicosities without significant edema, post sclerotherapy, small varicose dilations of cutaneous veins, initial varices during pregnancy.


Class 2: (23-32 mmHg)

Moderate symptoms of venous diseases such as aching heaviness and fatigue in the legs, chronic venous insufficiency, post sclerotherapy, prophylaxis of thrombosis, moderate to severe varicosities, post phlebitic syndrome, lymphedema, stasis dermatitis due to venous insufficiency, post fracture or post surgical edema, edema from burns, varices during pregnancy, pregnant patients with previous phlebitis.



AVAILABLE STYLES WITH

OPEN TOE



AD



AF*



AG*



AGG**



AT



ATM



ATU


	KL.1	KL.2	KL.1	KL.2	KL.1	KL.2	KL.1	KL.2	KL.1	KL.2	KL.1	KL.2
AD	162	262	AF*	162	262	AG*	162	262	AGG**	162	262	262
AT	162	262	ATM	162	262	ATU	162	262				

(34) Skin

* AF & AG available also with silicon band.

** AGG: adjustable waist belt, left or right style.

MEASURING POINTS





CIRCUMFERENCE MEASUREMENTS (cm)

	XS		S		M		L		XL	
	normal	max	normal	max	normal	max	normal	max	normal	max
g	42-46	47-54	47-51	52-60	52-57	58-67	58-63	64-74	64-70	71-83
f	37-40	41-47	41-45	46-52	46-50	51-59	51-56	57-65	57-62	63-73
d	28-30	31-36	31-34	35-40	35-38	39-44	39-42	43-50	43-47	48-55
b	18-19		20-22		23-25		26-28		29-31	

SIZES (cm)

XS - S - M - L - XL
NORMAL - MAX
SHORT - LONG
CUSTOM MADE

CUSTOM MADE (cm)



LENGTH MEASUREMENTS (cm)

	SHORT	LONG
A-D	34-38	39-44
A-F	52-59	60-68
A-G	62-71	72-82



Stay-up band with silicon drops



COMPOSITION: 66% Cotton - 19% Elastane - 15% Polyamide





LATEX FREE
CL. 2 & 3

UNIQUE STRETCH CHARACTERISTICS

The optimum elasticity allow a constant compression over long periods of time. The special porous mesh assures a normal transpiration of the skin. A complete range of styles and sizes, including Custom-made stockings, fit most patient needs.

COMPRESSION CLASSES:

- **Strong 261:** compression class 2
- **Strong 361:** compression class 3

INDICATIONS:

Class 2: (23-32 mmHg)

Moderate symptoms of venous diseases such as aching heaviness and fatigue in the legs, chronic venous insufficiency, post sclerotherapy, prophylaxis of thrombosis, moderate to severe varicosities, post phlebitic syndrome, lymphedema, stasis dermatitis due to venous insufficiency, post fracture or post surgical edema, edema from burns, varices during pregnancy, pregnant patients with previous phlebitis.


Class 3: (34-46 mmHg)

More severe presentations of indications for class 2; correctable lymphedema, severe chronic venous insufficiency, post phlebitic syndrome, chronic post surgical venous insufficiencies, severe post thrombotic syndrome, venous ulcers.



AVAILABLE STYLES WITH

OPEN TOE



	AD	AF*	AG*	AGG**	AT	ATM	ATU
KL.2	261	261	261	261	261	261	261
KL.3	361	361	361	361	361		361

(34) Skin


* AF & AG available also with silicon band.

** AGG: adjustable waist belt, left or right style.

SIZES (cm)

XS - S - M - L - XL
NORMAL - MAX
SHORT - LONG
CUSTOM MADE

MEASURING POINTS



CIRCUMFERENCE MEASUREMENTS (cm)


	XS		S		M		L		XL	
	normal	max	normal	max	normal	max	normal	max	normal	max
g	42-46	47-54	47-51	52-60	52-57	58-67	58-63	64-74	64-70	71-83
f	37-40	41-47	41-45	46-52	46-50	51-59	51-56	57-65	57-62	63-73
d	28-30	31-36	31-34	35-40	35-38	39-44	39-42	43-50	43-47	48-55
b	18-19		20-22		23-25		26-28		29-31	

LENGTH MEASUREMENTS (cm)


	SHORT	LONG
A-D	34-38	39-44
A-F	52-59	60-68
A-G	62-71	72-82

CUSTOM MADE (cm)

CIRCUMFERENCE MEASUREMENTS (cm)



LENGTH MEASUREMENTS (cm)



Stay-up band with silicon drops



COMPOSITION: Strong 261/361: 70% Polyamide - 30% Rubber Latex Free





COMPRESSION THERAPY AFTER LEG OPERATIONS

Kit indicated for postoperative surgery at the lower limbs consisting of an open toe thigh stocking with adjustable belt and a light understocking. The remarkable elasticity of the latex-free rubber and the special porous mesh assure effective therapeutic action and good skin transpiration throughout the day. The one size and light understocking helps to keep in place wound dressing facilitating the application of the elastic stocking.

Thigh stocking with increased compression at the thigh compared to a traditional compression stocking.

COMPRESSION CLASSES:

- Kit PostOp: compression classes 2

INDICATIONS:

For compression therapy after leg operations (sclerotherapy treatment, stripping, laser, radiofrequency, varicectomy).

PACKAGE CONTENTS:

- 1 thigh high with adjustable belt STRONG 261 (left or right)
- 2 understockings GLORIA FIX. Closed toe, one size.



AVAILABLE STYLES WITH
OPEN TOE



AGG**

CL.2


261

(34) Skin

** AGG style, thigh with adjustable belt left or right

The AG understockings (GLORIA FIX) to keep in place wound dressing are also available separately: single or pack with 20 pieces.

MEASURING POINTS



SIZES (cm)

XS - S - M - L - XL
NORMAL - MAX
SHORT - LONG
RIGHT - LEFT

CIRCUMFERENCE MEASUREMENTS (cm)

	XS		S		M		L		XL	
	normal	max	normal	max	normal	max	normal	max	normal	max
g	42-46	47-54	47-51	52-60	52-57	58-67	58-63	64-74	64-70	71-83
d	28-30	31-36	31-34	35-40	35-38	39-44	39-42	43-50	43-47	48-55
b	18-19		20-22		23-25		26-28		29-31	

LENGTH MEASUREMENTS (cm)

	SHORT	LONG	EXTRA-LONG
A-G	62-71	72-82	83-94

COMPOSITION: Monocollant Strong: 63% Polyamide - 37% Elastane
Understockings: 100% Polyamide



FLAT-KNIT MEDICAL COMPRESSION

The special texture of flat-knit garments makes them suitable for compression therapy in phlebology, lymphology and scar therapy. The unyielding compression ability of flat-knit, together with the massage effect produced, encourages lymphatic drainage and reduces edema expansion. The flat and elastic linked seams fit the shape of the body perfectly and make them suitable for patients with either deformed limbs, extra large proportions or medical complications such as elephantiasis.

CHARACTERISTICS:

- Soft flat-knit fabric
- Flat and elastic seams
- Comfortable and easy to wear
- 3 classes of compression (CL.2, CL.3, CL.4)
- Available in ready-to-wear standard sizes and made-to-measure
- Tailored foot caps

INDICATIONS:

- **Class 2:** Primary and secondary lymphedema (stage 1), lipoedema (stage 1), moderate post-surgical and post-traumatic edema, scar therapy. Surface inflammation of the veins after venous ulcers healing.
- **Class 3:** Primary and secondary lymphedema (stage 1 and 2), lipoedema (stage 1 and 2), scar therapy. Severe primary varicose veins with a tendency to edema and post-thrombotic syndrome.
- **Class 4:** Primary and secondary lymphedema (stage 2 and 3), lipoedema (stage 2 and 3), elephantiasis, serious vein disease with a tendency to edema, post-surgical and post-traumatic edema, or mixed forms such as phlebolymphe-
dema and lipo-lymphedema.

Absolute contraindications:

Infectious dermatitis, acute skin-hypodermatitis, dermatoses with ulcers or active sores. Arteriopathies with pressure index (ABI or IW) ≤ 0.55 (unless medically indicated), severe neuropathies, extrinsic vein compression (Baker's cyst of the popliteal cables, lymphadenopathy etc), decompensated heart failure, septic phlebitis, ischemic thrombosis (Phlegmasia Alba and Cerulea Dolens). Rheumatic fibromyalgia or generally any form of intolerance to prolonged pressure on the skin surface.


Relative contraindications:

Loss of limb sensitivity (unless medically indicated), intolerance to the material of compression stockings, advanced stage peripheral neuropathy (eg. Diabetes mellitus), chronic primary polyarthritis.



AVAILABLE STYLES STANDARD SIZES

OPEN TOE



AD

CL.2

CL.3

CL.4



AG

CL.2

CL.3



AGH

CL.2

CL.3



AGG*

CL.2

CL.3



AGT*

CL.2

CL.3


(34) Skin

(60) Black

*Models AGG and AGT, right or left

LEG CIRCUMFERENCE MEASUREMENTS (cm)

MEASURING POINTS



	4	5	6	7	8	9	10	11	12	13
g	44-46	46-48	48-50	50-52	52-54	54-56	56-58	58-60	60-62	62-64
f	41-42,5	42,5-44	44-45,5	45,5-47	47-48,5	48,5-50	50-51,5	51,5-53	53-54,5	54,5-56
e	31-32,5	32,5-34	34-35,5	35,5-37	37-38,5	38,5-40	40-41,5	41,5-43	43-44,5	44,5-46
d	28-29,5	29,5-31	31-32,5	32,5-34	34-35,5	35,5-37	37-38,5	38,5-40	40-41,5	41,5-43
c	30-31,5	31,5-33	33-34,5	34,5-36	36-37,5	37,5-39	39-40,5	40,5-42	42-43,5	43,5-45
b ¹	24-25,5	25,5-26,5	26,5-27,5	27,5-29	29-30	30-31,5	31,5-32,5	32,5-34	34-35	35-36
b	19-20	20-21	21-22	22-23	23-24	24-25	25-26	26-27	27-28	28-29
h	28-29	29-30	30-31	31-32	32-33	33-34	34-35	35-36	36-37	37-38
a	19-20	20-21	21-22	22-23	23-24	24-25	25-26	26-27	27-28	28-29

LEG LENGTH MEASUREMENTS (cm)

/D = 39-41

/G = 74-78

COMPOSITION: Polyamide - Elastane



AVAILABLE STYLES MADE TO MEASURE

OPEN TOE



AD

CL.2

CL.3

CL.4



AG

CL.2

CL.3

CL.4



AGH

CL.2

CL.3

CL.4



AGG*

CL.2

CL.3

CL.4



AGT*

CL.2

CL.3

CL.4



AT

CL.2

CL.3

CL.4



ATM

CL.2

CL.3

CL.4



ATU

CL.2

CL.3

CL.4

(34) Skin


(60) Black

*Models AGG and AGT, right or left


VARIATIONS FOR MADE TO MEASURE




diagonal toe



diagonal thigh



knee bend



closed toe

VARIATIONS FOR MADE TO MEASURE LEG

- Complete compression
- Higher compression above
- Lower compression below.....

VARIATIONS FOR MADE TO MEASURE TIGHTS

- Complete compression
- Different compressions (specify)
- ATM with front flap or open

FOOT CAP MADE TO MEASURE

AVAILABLE STYLES



open toe - 4 toes

CL.2

CL.3



open toe - 5 toes

CL.2

CL.3



closed toe - 4 toes

CL.2

CL.3



closed toe - 5 toes

CL.2

CL.3

(34) Skin

(60) Black

MEDICAL COMPRESSION STOCKINGS IN SYNTHETIC, SOFT, LIGHTWEIGHT AND EASY TO WEAR

Easyven® stockings, in soft and elastic synthetic yarns, combine fit and comfort without losing the therapeutic efficacy. Easy to put on and remove provide comfort and well-being throughout the day.

COMPRESSION CLASSES:

- EasyVen: 20-30 mmHg


INDICATION:

moderate Chronic venous insufficiency (CVI) and after: healing of superficial ulcer, phlebitis syndrome, sclerotherapy treatment, ambulant selective varicose surgery, stripping. Varices during pregnancy, prophylaxis of thrombosis/embolism.

CONTRAINDICATIONS:

infections, eczema, ulcers of lower limbs of arterial origins, severe neuropathies, uncompensated heart failure, severe leg edema, pulmonary edema, pre-existing gangrene, severe arterial circulation disorders.

Stay-up band with silicon drops



*AG available also with stay-up silicon band

**AGG: adjustable waist belt, left or right style

AVAILABLE STYLES:

OPEN TOE



AD




AG*



AGG**

(34) Skin

MEASURING POINTS



CIRCUMFERENCE MEASUREMENTS (cm)

	S		M		L		XL	
	normal	max	normal	max	normal	max	normal	max
g	47-51	52-60	52-57	58-67	58-63	64-74	64-70	71-83
d	31-34	35-40	35-38	39-44	39-42	43-50	43-47	48-55
b	20-22		23-25		26-28		29-31	

LENGTH MEASUREMENTS (cm)

	SHORT	LONG
A-D	34-38	39-44
A-G	62-71	72-82

COMPOSITION: Polyamide - Elastane



ULCERKIT®**GLORIA ULCERKIT®
THE REVOLUTION IN THE
TREATMENT OF VENOUS ULCERS.**

Ulcerkit® is a therapeutic compression stocking system indicated for leg ulcers, easy to use improves patients' compliance. Gloria Ulcerkit® can also be used as a normal medical stocking when a certain effort for application is required.

COMPRESSION CLASSES:

UlcerKit® is available in Class 1 and 2. As the system is composed by two separated parts, in case of need of a higher compression near the ankle (where usually the ulcer is located), there is the possibility to combine the lower part in KKL. 2 with the upper part in KKL. 1.

So the range is the following:

- **UlcerKit® 20/20:**
ankle KKL. 1 - calf KKL. 1
- **UlcerKit® 30/20:**
ankle KKL. 2 - calf KKL. 1
- **UlcerKit® 30/30:**
ankle KKL. 2 - calf KKL. 2

Patented Product IT 1319230 – EU 1337211



INDICATIONS:**Class 1:**

treatment of venous ulcer associated with other compression aids. Chronic venous insufficiency, post phlebitic, post ulcer, after sclerotherapy. Prevention of recidivism. Initial varices during pregnancy.

Class 1/2:


treatment of venous ulcer when a higher ankle compression is required. Stronger action at the ankle and moderate at the calf.

Classe 2:

same as class 1 and 1-2, but with a highest compression on the leg, edema.

CONTRA-INDICATIONS:

Do not wear UlcerKit® in presence of severe arterial venous insufficiency. In moderate cases UlcerKit® can be used and in case of pain the stocking can be removed.

AVAILABLE STYLE: AD**CIRCUMFERENCE MEASUREMENTS (cm)**

	S		M		L		XL	
	normal	max	normal	max	normal	max	normal	max
c	30-33	34-37	34-37	38-41	38-41	42-45	42-45	46-49
b	19-21		22-24		25-27		28-31	

POINT B**POINT C****HOW TO WEAR THE KIT:**

Step 1



Step 2



Step 3



Step 4



Step 5





GLORIA ULCERKIT® PRO MULTI-LAYER STOCKING SYSTEM

Gloria UlcerKit PRO® is a multi-layer stocking system to treat venous ulcers and is available in two different compressions:

- UlcerKit PRO 35 STRONG (35-40 mmHg)
- UlcerKit PRO 35 COTTON (35-40 mmHg)
- UlcerKit PRO 40 STRONG (40-45 mmHg)
- UlcerKit PRO 40 COTTON (40-45 mmHg)


INDICATIONS:

UlcerKit® PRO 35 (35–40 mmHg):
skin ulcers of the leg without venous insufficiency (eg post-traumatic, vasculitic), skin ulcers of the leg with mild to moderate venous insufficiency.

UlcerKit® PRO 40 (40–45 mmHg):
skin ulcers of the leg with severe venous insufficiency.

CONTRAINDICATIONS:

Gloria UlcerKit PRO® should not be used in the presence of significant arterial disease with ankle-brachial index (ABI) of less than 0.8. Diabetes mellitus with diabetic neuropathy.



PACKAGE CONTENTS:

• UlcerKit® PRO 35

Calf-length KL.2* (1 pcs)
Calf-length 18 mmHg (1 pcs)
Thin understocking AD (2 pcs)
Foot slip


• UlcerKit® PRO 40

Calf-length KL.2* (1 pcs)
Calf-length 24 mmHg (1 pcs)
Thin understocking AD (2 pcs)
Foot slip



* Cotton or Strong series

AVAILABLE STYLE: AD



CIRCUMFERENCE MEASUREMENTS (cm)

	S		M		L		XL	
	normal	max	normal	max	normal	max	normal	max
d	31-34	35-38	34-37	39-42	39-42	43-47	43-47	48-52
c	32-35	36-39	36-39	40-43	40-43	44-48	44-48	49-53
b	20-22		23-25		26-28		29-31	

LENGTH MEASUREMENTS (cm)

	SHORT	LONG
A-D	34-36	39-44

POINT B



POINT C



POINT D



HOW TO WEAR THE KIT:



Step 1



Step 2



Step 3



Step 4



Step 5



Step 6



Step 7



18 mmHg



THROMBOEMBOLIC PROPHYLAXIS

A complete line of anti-embolism stockings. Comfortable and easy to wear thanks to a very soft fabric, well tolerated due to a normal transpiration of the skin. Compression of 18 mmHg at the ankle.

COMPRESSION CLASS:

- **THROMBOSIS 18 SYNTHETIC:** 18 mmHg
- **THROMBOSIS 18 COMFORT COTTON:** 18 mmHg

AVAILABLE STYLES WITH:

CLOSED TOE with INSPECTION HOLE



AD

18 mmHg



AG*

18 mmHg



AGG

left or right

18 mmHg



AGG**

ambidextrous

18 mmHg

(10) White

* AG available also with silicon band.


** AGG available only for Thrombosis 18 Synthetic

INDICATIONS

Clinically designed to promote venous return for prevention of thrombosis before, during and after surgery and for medical conditions requiring prolonged bed rest. Before, during and after childbirth. Bedridden oncologic patients.

Contraindications: Infections, eczema, ulcers of lower limbs, severe neuropathies, uncompensated heart failure, severe leg edema, pulmonary edema, pre-existing gangrene, severe arterial circulation disorders.

MEASURING POINTS



CIRCUMFERENCE MEASUREMENTS (cm)

	S		M		L		XL		XXL
	normal	max	normal	max	normal	max	normal	max	normal
g	45-54	53-62	53-62	61-70	61-70	69-78	69-78	77-86	77-86
c	30-36		35-41		40-46		45-51		50-56
b	19-21		22-24		25-27		28-30		31-33

LENGTH MEASUREMENTS (cm)

	SHORT	REGULAR	LONG
A-D	35-39	39-44	44-50
A-G	63-72	72-82	82-93

COMPOSITION:

Thrombosis 18:

Polyamide - Elastane

Thrombosis 18 Comfort:

Polyamide - Cotton - Elastane

POSTOP®

24 mmHg

POST-SURGERY STOCKINGS


For Venous thromboembolism (VTE) prophylaxis in the post operative phase in patients at risk or with a history of previous thrombosis. Providing 24 mmHg compression, Gloria Post-Op series apply a higher pressure than the traditional anti-thrombosis stockings. Useful in vascular surgery, orthopedic and cardiac surgery.

COMPRESSION CLASS:

- 24 mmHg

AVAILABLE STYLES WITH:

OPEN TOE



AD

24 mmHg



AG*

24 mmHg



AGG**

24 mmHg

(10) White

* AG available also with silicon band.


** AGG: adjustable waist belt, left or right style.

INDICATIONS

For partially mobilized patients before and after surgery, for post venous stripping, for those who have already had a thrombosis, and for immobilized patients with venous diseases. Postoperative edema in orthopedic surgery.

Contraindications: Infections, eczema, ulcers of lower limbs, severe neuropathies, uncompensated heart failure, severe leg edema, pulmonary edema, pre-existing gangrene, severe arterial circulation disorders.

MEASURING POINTS



CIRCUMFERENCE MEASUREMENTS (cm)

	S		M		L		XL		XXL
	normal	max	normal	max	normal	max	normal	max	normal
g	45-54	53-62	53-62	61-70	61-70	69-78	69-78	77-86	77-86
c	30-36		35-41		40-46		45-51		50-56
b	19-21		22-24		25-27		28-30		31-33

LENGTH MEASUREMENTS (cm)

	SHORT	LONG
A-D	34-39	39-45
A-G	62-72	72-84

COMPOSITION:

Post-Op 24:
Polyamide - Elastane





GLORIA MED® SLEEVE 101 – 201

Elastic compression arm-sleeves GloriaMed® lightweight and comfortable, thanks to the use of fine synthetic yarns, provide excellent comfort to the superior limb.

COMPRESSION CLASSES:

Gloria Med® Sleeve 101:

Compression classes 1

Gloria Med® Sleeve 201:

Compression classes 2

FEATURES

- Available styles: with or without mitten
- Standard sizes and custom-made
- 2 Lengths (short or long)
- Soft and comfortable to wear
- Thin and flattened seams

INDICATIONS:

Arm-sleeve are indicated in the treatment of primary and secondary lymphedema of the upper limbs.

Class 1:

prevention of post-surgery lymphedema after mastectomy with or without chemo and radiological treatments. Treatment of mild lymphedema. Prophylaxis of recurrent edema after decongestive therapy.

Class 2:

Prevention of edema after decongestion is achieved. Post-mastectomy fluid accumulation and edema following radiation therapy. Post-surgery and post-traumatic edema after basic treatment.

CONTRAINDICATIONS:


- Uncontrolled lymphedema
- Inflammation
- Open wounds
- Infections
- Arterial circulatory problems of the upper limb




AVAILABLE STYLES

WITH SHOULDER AND STRAP

WITH OR WITHOUT SILICON BAND



CG



AG



CG



AG

KL.1

101

101

101

101

KL.2

201


201

201

201

(34) Skin

MEASURING POINTS



SIZES (cm)

S - M - L - XL
SHORT - LONG
CUSTOM MADE

CIRCUMFERENCE MEASUREMENTS (CM)


	S	M	L	XL
g	25-30	29-34	33-38	37-42
d	21-25	24-28	27-31	30-34
c	15-17	17-19	19-21	21-23

LENGTH MEASUREMENTS


	SHORT	LONG
C-G	41-46	46-52

CUSTOM MADE (cm)

Circumference measurements



Length measurements





FLAT KNIT SLEEVE AND GLOVE

Ideal for upper limb lymphedema



The special texture of flat-knit garments makes them suitable for compression therapy in phlebology, lymphology and scar therapy. The unyielding compression ability of flat-knit, together with the massage effect produced, encourages lymphatic drainage and reduce edema expansion. The flat and elastic seams fit the shape of the body perfectly.

CHARACTERISTICS:

- Soft flat-knit fibre
- Flat and elastic seams
- Comfortable and easy to wear
- 2 classes of compression (CL. 2, CL. 3)
- Sleeves and gloves available in standard sizes and made to measure with different variations available

INDICATIONS

Primary and secondary lymphedema of the upper limb (stage 2 and 3), lipoedema (stage 2 and 3), post surgical and post traumatic edema, upper extremity deep vein thrombosis (Paget von Schrötter syndrome), hypertrophic scars (following burns for example), rheumatoid arthritis.

Absolute contraindications:

Ischemia (eg. Severe arterial insufficiency), decompensated heart failure, untreated septic phlebitis, ischemic thrombosis (Phlegmasia coerulea dolens). Rheumatic fibromyalgia or generally any form of intolerance due to prolonged pressure on the skin surface. Skin infections, acute or exudative dermatitis.

Relative contraindications:

Loss of limb sensitivity (unless medically indicated), intolerance to compression stocking material.


ARM CIRCUMFERENCE MEASUREMENTS (cm)

	1		2		3		4	
	normal	max	normal	max	normal	max	normal	max
g	< 32	< 36	32 - 36	36 - 40	36 - 40	40 - 44	> 40	> 44
f	< 27	< 29	27 - 30	29 - 32	30 - 33	32 - 35	> 33	> 35
e	< 25		25 - 28		28 - 30		> 30	
d	< 24		24 - 26		26 - 29		> 29	
c	< 18		18 - 20		20 - 22		> 22	
b	< 19		19 - 21		21 - 23		> 24	
a	< 18		18 - 20		20 - 22		> 22	

LENGTH MEASUREMENTS (cm)






	SHORT	LONG
c-g	< 44	> 44

MEASURING POINTS



FLAT KNIT SLEEVE **STANDARD** AND **MADE TO MEASURE**






AVAILABLE MODELS

				
CG without mitten	AG with mitten	AG mitten with thumb	AG with open finger tips	AG with closed finger tips
CL.2	CL.2	CL.2	CL.2	CL.2
CL.3	CL.3	CL.3	CL.3	CL.3

(34) Skin





(60) Black

SLEEVE VARIATIONS

				
3 cm silicon border	5 cm silicon border	Diagonal top border	Shoulder cap with adjustable strap	Shoulder cap with velcro bra attachment

STANDARD**MADE TO MEASURE**FLAT KNIT GLOVES **STANDARD** AND **MADE TO MEASURE**

AVAILABLE MODELS

			
mitten	mitten with thumb	mitten with open finger tips	mitten with closed finger tips
CL.2	CL.2	CL.2	CL.2
CL.3	CL.3	CL.3	CL.3

(34) Skin

(60) Black



ELEGANCE AND MEDICAL COMPRESSION IN FASHIONABLE COLOURS

Attractive, sheer and comfortable stockings and pantyhose with graduated compression. Suitable for young and dynamic women who need a medical effect with elegance and style.

INDICATIONS

8-10 mmHg:

Support pantyhose physician-recommended for prolonged standing, sitting and for young ladies during pregnancy.

12-15 mmHg:

Mild compression for tired and heavy legs. Prevention and support in presence of mild varicosity during pregnancy.

18-20 mmHg:

Designed specially to provide therapeutic benefits and compression in presence of heaviness and fatigue in the legs, initial varices during pregnancy.

20-30 mmHg:

Mild to moderate symptoms of venous diseases such as heaviness and fatigue in the legs, post-sclerotherapy and after surgery, mild to moderate varicosities without significant edema, in presence of more severe varices during pregnancy.




LADYGLORIA 8

8-10 mmHg - SHEER - SIZE 1-5

AVAILABLE STYLES WITH

CLOSE TOE

OPEN TOE



AD



AG


AG
grip top


AT



ATM



AD


AG
grip top


AT

White

Sahara

Deer

Sierra

Beaver

Blu Marine

Dark grey

Black

LADYGLORIA 12

12-15 mmHg - FISH NET - SIZE 1-5

AVAILABLE STYLES WITH

CLOSE TOE

OPEN TOE



AD



AG


AG
grip top


AT



ATM



AD


AG
grip top


AT

White

Sahara

Deer

Sierra

Beaver

Safari

Blu Marine

Dark grey

Black




LADYGLORIA 18

18-20 mmHg - FISH NET - SIZE 1-5

AVAILABLE STYLES WITH

CLOSE TOE

OPEN TOE



AD



AG




AG
grip top



AT



ATM



AD



AG
grip top



AT

WHITE

Sahara

Deer

Sierra

Beaver

Blu Marine

Dark grey

Black



LADYGLORIA 24

20-30 mmHg - FISH NET - SIZE 1-6

AVAILABLE STYLES WITH

CLOSE TOE

OPEN TOE



AD




AG
grip top



AT



ATM



AD



AG
grip top



AT

Deer

Black



SIZE TABLE FOR: AG, AT, ATM

		Weight Kg																			
		45	48	50	52	54	57	59	61	64	66	68	70	73	75	77	79	81	84	86	88
Height cm	150																				
	152																				
	155																				
	157																				
	160																				
	163																				
	165																				
	168																				
	170																				
	173																				
	175																				
	178																				
	180																				
	182																				

TABLE FOR: AD

SIZES	SHOE SIZE	
	EU	USA
1	35 - 36	5 - 6
2	37 - 38	6½ - 7½
3	38,5 - 39	8 - 8½
4	39,5 - 40	9 - 9½
5	41 - 42	9½ - 10

SIZES (cm)

1 - 5

*Lady Gloria 24: 1 - 6

COMPOSITION: Polyamide - Elastane



A COMPLETE RANGE FOR DYNAMIC MEN

Preventive socks of high quality, available in several colours and 2 compression selections. Their elegance and comfort permit wearing them throughout the day whether at work or leisure.

The cotton style provides increased transpiration and reduced allergies.

INDICATIONS:

15-18 mmHg:

Support socks physician-recommended for prolonged standing, sitting or travelling. For tired and heavy legs.


18-20 mmHg:

Designed specially to provide therapeutic benefits and compression in presence of heaviness and fatigue in the legs. Prophylaxis of venous disorders in patients at risk.



MAN'SUPPORT 15

15-18 mmHg




AD

	SINTHETIC	COTTON
White	<input type="checkbox"/>	<input type="checkbox"/>
Black	<input type="checkbox"/>	<input type="checkbox"/>
Dark Grey	<input type="checkbox"/>	<input type="checkbox"/>
Brown	<input type="checkbox"/>	<input type="checkbox"/>
Bordeaux	<input type="checkbox"/>	<input type="checkbox"/>
Blue	<input type="checkbox"/>	<input type="checkbox"/>

MAN'SUPPORT 18

18-20 mmHg



AD

	SINTHETIC	COTTON
White	<input type="checkbox"/>	<input type="checkbox"/>
Black	<input type="checkbox"/>	<input type="checkbox"/>
Dark Grey	<input type="checkbox"/>	<input type="checkbox"/>
Brown	<input type="checkbox"/>	<input type="checkbox"/>
Bordeaux	<input type="checkbox"/>	<input type="checkbox"/>
Blue	<input type="checkbox"/>	<input type="checkbox"/>

SIZES	SHOE SIZE	
	EU	USA
1	36 - 38	5 - 6
2	39 - 40	6½ - 7
3	41 - 42	8 - 8½
4	43 - 44	9½ - 10
5	45 - 46	11 - 11½
6*	47 - 48	13 - 14

* Available Size only for ManSupport 18 Cotton blue and black.

COMPOSITION: Man'Support: Polyamide - Elastane
Man'Support Cotton: Polyamide - Cotton - Elastane



MAN'SUPPORT[®] 18

SUPPORT SOCKS

SILVER COTTON




SUPPORT SOCKS FOR MEN WITH SILVER IONS

Attractive and comfortable cotton socks with silver ions. Thanks to antibacterial, odor-control, thermoregulation properties of silver they are suitable for all occasions. Physician-recommended for tired and heavy legs.

INDICATIONS:


Designed specially to provide therapeutic benefits and compression in presence of heaviness and fatigue in the legs. Prophylaxis of venous disorders.



MAN'SUPPORT 18 — 18-20 mmHg



AD



SIZES	SHOE SIZE	
	EU	USA
1	36 - 38	5 - 6
2	39 - 40	6½ - 7
3	41 - 42	8 - 8½
4	43 - 44	9½ - 10
5	45 - 46	11 - 11½
6	47 - 48	13 - 14

COMPOSITION: Polyamide (Q-Skin[®]) - Cotton - Elastane



DIABE Soft®



COMFORT AND PROTECTION FOR DIABETIC AND SENSITIVE FEET

DIABESoft® is a line of socks indicated for people who suffer from diabetes, Athlete's Foot, sensitive and delicate feet. The special cotton yarn with integrated silver has an antibacterial effect, anti-odor and regulate the temperature keeping the foot healthy, dry and protected. Without compression to not compress and avoid tourniquet effect. Without seams and with soft terry insole, to avoid skin abrasions and ensure a perfect fit and maximum comfort.



DIABESoft

AVAILABLE STYLES:

CLOSED TOE



A-B1



AD

Black

Skin

SIZES	SHOE SIZE
S	35 - 37
M	38 - 40
L	41 - 43
XL	44 - 46
XXL	47 - 48

COMPOSITION: Cotton (with integrated silver) - Elastane





GRADUATED COMPRESSION SOCKS FOR MAN & WOMAN

Sitting or standing for long periods of time without moving your legs is a condition that causes a constriction of the veins and a slow return of venous blood to the heart. Swollen and painful legs are the first possible signs.

The graduated compression socks for men and women "Travel & Work" contribute actively to the health of your legs.

INDICATIONS:

15-20 mmHg:

Ideal for long trips by car, plane, train, and for those who in everyday life maintains a static position for many hours, or with their legs crossed.



TRAVEL & WORK

15-20 mmHg



AD

Black

Blue

TABLE FOR: AD

SIZES	SHOE SIZE	
	EU	USA
1	36 - 38	5 - 6
2	39 - 40	6½ - 7
3	41 - 42	8 - 8½
4	43 - 44	9½ - 10
5	45 - 46	11 - 11½

COMPOSITION: Polyamide - Elastane





FRESHNESS, COMFORT AND RELIEF OF TIRED LEGS

A Derma cosmetic formulation made by functional ingredients of natural origin, like Melilot, Centella, Rice bran oil, Alpha-Bisabolol, Vitamin E, Menthyl Lactate, that give to your legs' skin an immediate feeling of freshness, comfort and relief.

INDICATIONS:

- Topical adjuvant in the treatment of limphedema and post-traumatic edema.
- Venous insufficiency
- Swelling, heaviness, fatigue, cramps, sensation of warmth or burning, itching and tingling of the legs.
- To improve superficial blemish and improve microcirculation
- Indicated after sclerotherapy and laser therapy.

Tube 150 ml

BIBLIOGRAPHY:

- Judit R. Casley-Smith, John R. Casley-Smith; Modern treatment of lymphoedema. II. The benzopyrones. *Australas J Dermatol.* 1992; 33(2): 69-74
- Judit R. Casley-Smith, John R. Casley-Smith; Treatment of lymphoedema by complex physical therapy, with and without oral and topical benzopyrones: what should therapists and patients expect. *Lymphology*, 1996 Jun;29 (2): 76-82
- Cesarone MR, Laurora G, De Sanctis MT, Belcaro G.; Activity of Centella asiatica in venous insufficiency. *Minerva Cardioangiol.* 1992 Apr; 40(4):137-43
- Bell DR, Gochenaur K; Direct vasoactive and vasoprotective properties of anthocyanin-rich extracts. *Appl Physiol*; 2006 Apr; 100(4):1164-70. Epub 2005 Dec 8.



VITALITY AND SOFTNESS FOR YOUR LEGS

A Derma cosmetic formulation made by functional ingredients of natural origin, as Centella, Black Blueberry, Vitamin PP, Rice bran oil, Alpha-Bisabolol, Vitamin E, Menthyl Lactate, that give to dried skin of the legs an immediate vitality and softness.

INDICATIONS:

- Venous Insufficiency
- Especially for dry skin of the legs, thin and delicate
- Swelling, heaviness, fatigue.
- To improve superficial blemish and improve microcirculation
- Indicated after Sclerotherapy and laser therapy in subjects with sensitive skin.

Tube 150 ml

BIBLIOGRAPHY:

- Cesarone MR, Laurora G, De Sanctis MT, Belcaro G.; Activity of Centella asiatica in venous insufficiency. *Minerva Cardioangiol.* 1992 Apr; 40(4):137-43
- Bell DR, Gochenaur K; Direct vasoactive and vasoprotective properties of anthocyanin-rich extracts. *Appl Physiol*; 2006 Apr; 100(4):1164-70. Epub 2005 Dec 8.



USE: We recommend 2 applications a day, preferably one in the morning and one in the evening. It can be applied before to wear medical stockings as it does not grease, doesn't stain and absorbs quickly.

WARNINGS: Only for external use. Avoid contact with eyes and mucous membranes. Do not apply on wounds and injuries. Keep out of reach of children

VALIDITY: 36 months if not opened. Validity post-opening 12 months (PAO).





GEL FOR THE TREATMENT OF SKIN PIGMENTATIONS OF HEMOSIDERIN ORIGINS

Siderlink is a patented thioglycolic acid-based formulation at 20% (2-MERCAPTOCETIC ACID AT 20%) The product efficacy is due to its capacity to solubilize the hemosiderin and to bond the iron ion using a single active ingredient (2-mercaptoacetic acid).

INDICATIONS:

Treatment of cutaneous hyperpigmentations of hemosiderin origins which can often occur during CVI (chronic venous insufficiency), in Phlebology after sclerotherapy, laser, varicectomy. After post-traumatic and surgical therapies or aesthetical treatments (dermo-abrasion, peeling etc). Intravenous infusion administrations.

The product is also useful on epidermal spots of melanic origin.

THE PACKAGE CONTAINS:

- **SIDERLINK** gel (thioglycolic acid at 20%)
15 single use tubes 2 g/each
- **DETERLINK** Cleansing milk to clean the skin
1 vial of 50 ml
- **NEUTROLINK** Neutralizing lotion
2 vials of 50 ml/vial (Buffer solution to wash and neutralize the treated site)

WARNINGS: The use on the skin is generally well tolerated; but it is very important to test any skin reactivity during the first applications leaving the gel on the treated part for a short application time. Do not use it if there is a skin intolerance. Avoid eyes and mucous membrane contact and if it happens wash immediately the site with tap water or with weakly basic solutions such as bicarbonate solution at 5%. Do not ingest the product and keep it away from children.

CONTRAINDICATIONS: Do not use the product on skin lesions (abrasions, excoriations, burns, phlyctena), after sun exposition in summer time, just after wax depilatory treatments.

VALIDITY: 24 months for unopened package.

SIDERLINK (small tube): after opening use within 30 days.

DETERLINK/NEUTROLINK: PAO 12 months.

Once open keep the kit in a dry place, carefully close the cup of Siderlink after use to avoid the contact with air (the product is subject to easy oxidation).

CATEGORY OF M.D.: Siderlink is a medical device in class 1

Gloria Med Pharma S.r.l. patent EU 0975321



GLORIA MED[®] CUSTOM-MADE ORDER FORM

Fax or email to:



+39 0344 32035



gloriamed@gloriamed.it

Gloria Med S.p.A.

Custom Made Department

Via Diaz, 7 - 22017 Menaggio (Co)


Italy



= available



= not available



Compression	Item GLORIA-MED	Material	Toe	A-D	A-F	A-G	A-GH silicon band	A-GG	A-T	A-TM	A-TU
Class 1	COTTON 162	cotton	open					left right			
Class 2	STRONG 261	rubber	open								
	COTTON 262	cotton	open								
Class 3	STRONG 361	rubber	open								

Important: Indicate the exact circumference and length of the limb at each specific point. The measurements must be taken with limb resting and relaxed. For A-D, A-F and AG/AGH, please specify if a single item or pair is required.

☐ piece ☐ pair

CIRCUMFERENCE MEASUREMENTS (cm)		LENGTH MEASUREMENTS (cm)	
Waist	cT	IK ₁ T	IT
Hips	cH	IK ₁	IH
Left leg	cK	IG	IF
	cG	IE	ID
	cF	IC	IB ₁
	cE	IB	IA
	cD		
	cC		
	cB ¹		
	cB		
	cA		
Right leg	cK	IK ₂ T	IT
	cG	IK ₂	IH
	cF	IG	IF
	cE	IE	ID
	cD	IC	IB ₁
	cC	IB	IA
	cB ¹		
	cB		
	cA		

Type or write in block letters

DATE _____

PATIENT _____

PHYSICIAN _____

RETAILER _____

Stamp & Signature _____

Please write clearly

DATE: _____

PATIENT: _____



DOCTOR: _____




CLINIC _____

Address _____

Customer signature: _____

Note/Special requests/ Colour: _____

Choose article and compression	Quantity			Colour
FLAT KNIT	RIGHT	LEFT	TIGHTS	Skin <input type="checkbox"/> Black <input type="checkbox"/>
Gloriamed Perfect CCL. 2 <input type="checkbox"/>				
Gloriamed Perfect CCL. 3 <input type="checkbox"/>				
Gloriamed Perfect CCL. 4 <input type="checkbox"/>				
Toe seam				
OPEN TOE <input type="checkbox"/>				
CLOSED TOE <input type="checkbox"/>				
	Model THIGH with waist attachment  Closure: <input type="checkbox"/> Velcro <input type="checkbox"/> Adjustable buttons Waist circumference _____ Length G to waist (taken on hip) right _____ left _____			ONE LEGGED PANTY 

CIRCUMFERENCE	LENGTH	VARIATIONS
<div style="display: flex; justify-content: space-around;"> <div style="text-align: center;"> RIGHT LEG  </div> <div style="text-align: center;"> LEFT LEG  </div> </div> 	<div style="display: flex; justify-content: space-between;"> <div> Front PANTY SIZE _____ cm Pubic bone to t _____ cm g-t _____ (take measurements on hip) j-t _____ (take measurements on hip) Right leg _____ A-g _____ _____ A-f _____ _____ A-e _____ _____ A-d _____ _____ A-c _____ _____ A-b1 _____ _____ A-b _____ Right foot Foot length (h-a) _____ cm Inside foot length _____ cm (for closed toe models) Left foot Foot length (h-a) _____ cm Inside foot length _____ cm (for closed toe models) </div> <div> Back _____ cm Gluteal fold to t _____ Left leg _____ A-g _____ _____ A-f _____ _____ A-e _____ _____ A-d _____ _____ A-c _____ _____ A-b1 _____ _____ A-b _____ </div> </div>	PANTY VARIATIONS Full compression <input type="checkbox"/> Different compression <input type="checkbox"/> (indicate in above notes) Open front for men <input type="checkbox"/> Without gusset <input type="checkbox"/> Maternity: with front <input type="checkbox"/> open front <input type="checkbox"/> LEG VARIATIONS Full compression <input type="checkbox"/> Higher compression above <input type="checkbox"/> Reduced compression above <input type="checkbox"/> Diagonal thigh (length A-g1) <input type="checkbox"/> Knee bend <input type="checkbox"/> LACE HOLD-UP 5 cm of silicon dots <input type="checkbox"/> FOOT OPENING Diagonal foot (open toe) <input type="checkbox"/> Internal length foot _____ cm External length of foot _____ cm

Please write clearly

DATE: _____

PATIENT: _____

DOCTOR: _____

CLINIC _____

Address _____

Customer signature: _____

Note/Special requests/ Colour: _____

Choose article:


GloriaMed Perfect CCL2 ☐

GloriaMed Perfect CCL3 ☐

CIRCUMFERENCE	TOES	Toe 5	Toe 4	Toe 3	Toe 2	Toe 1	TOES	Toe 1	Toe 2	Toe 3	Toe 4	Toe 5	CIRCUMFERENCE
	Circumf. base nail						Circumf. base nail						
	Circumf. base toe						Circumf. base toe						
	Toe length						Toe length						

LEFT FOOT

WITHOUT TOE 5 ☐




QUANTITY COLOUR

OPEN TOE ☐ CLOSED TOE ☐

LEFT FOOT	
LENGHT 1	cm
LENGHT 2	cm

RIGHT FOOT

WITHOUT TOE 5 ☐



QUANTITY COLOUR

OPEN TOE ☐ CLOSED TOE ☐

RIGHT FOOT	
LENGHT 1	cm
LENGHT 2	cm

Fax or email to:



+39 0344 32035







gloriamed@gloriamed.it

Gloria Med S.p.A.

Custom Made Department

Via Diaz, 7 - 22017 Menaggio (Co)

Italy

		with shoulder and strap		with silicon band	
					
Compression	Item	with mitten	without mitten	with mitten	without mitten
Class 1	ARM SLEEVE 101				
Class 2	ARM SLEEVE 201				

When a standard size is not possible a made to measure arm-sleeve is required. Please be so kind to fill in all arm circumferences and lengths.

☐ Right _____ ☐ pcs _____ ☐ Left _____ ☐ pcs _____

Circumference measurements

RIGHT

G-H cm

Belt Circumference cm

G-H cm

g cm

f cm

e cm

d cm

c1 cm

c cm

b cm

a cm

A-B cm

A-C cm

C-D cm

C-E cm

C-F cm

C-G cm

LEFT

G-H cm

g cm

f cm

e cm

d cm

c1 cm

c cm

b cm

a cm

A-B cm

A-C cm

C-D cm

C-E cm

C-F cm

C-G cm

Length measurements

ICG cm

ICF cm

ICE cm

ICD cm

ICC1 cm

IAC cm

IAB cm

A-B cm

A-C cm

C-D cm

C-E cm

C-F cm

C-G cm

Type or write in block letters

DATE _____ PATIENT _____

PHYSICIAN _____

RETAILER _____

Stamp & Signature

Please write clearly

DATE: _____

PATIENT: _____

DOCTOR: _____





CLINIC _____


Address _____








Customer signature: _____


Note/Special requests/ Colour: _____

☐ Black ☐ Skin

NUMBER OF PIECES	LEFT	RIGHT
Compression class	Glove model	
FLAT KNIT	 Hand mitten (no fingers with thumb hole)	
	 Hand mitten with thumb (no fingers)	
	 Hand mitten with fingers, open tips	
	 Hand mitten with fingers, closed tips	
Glove CCL2	<input type="checkbox"/>	
Glove CCL3	<input type="checkbox"/>	



NUMBER OF PIECES	LEFT	RIGHT
Compression class	Sleeve model	
FLAT KNIT	 Sleeve to wrist (c-g)	
	 Sleeve with mitten (a-g), with thumb hole and no fingers	
	 Sleeve with mitten (a-g), with thumb, no fingers	
	 Sleeve with mitten (a-g), with thumb and fingers	
Sleeve CCL2	<input type="checkbox"/>	
Sleeve CCL3	<input type="checkbox"/>	
Silicon band variations		
		
		
Silicon band 3 cm	<input type="checkbox"/>	
Silicon band 5 cm	<input type="checkbox"/>	
Sleeve variations	Finger type Open <input type="checkbox"/> Closed <input type="checkbox"/>	




Sleeve with diagonal top border

G - G' _____ cm



Sleeve with shoulder cap and velcro that attaches to bra

G - H _____ cm



Sleeve with shoulder cap and adjustable strap

G - H _____ cm



Gardamed Limited, The Pixnmore Centre,
Pixmore Avenue, Letchworth Garden City. Hertfordshire.
SG6 1JG. Tel. 01462411819 info@gardamed.com